



Instructions for Use

【Brand Name】 FriCARE

【Model No】 Ecoguard A

【Product Name】 Disposable Face Mask

【Procedure for Use】 Hold the color side of the face mask facing outwards, and place it over the nose and mouth. Secure ear loop around each ear. Mold the adjustable nose strip into the shape of your nose, and pull the bottom of the face mask under the chin. Discard after use.

【Intended Use】 The materials used for our disposable face masks meet the criteria of ASTM Level 3 performance requirements, and suitable for use by the general public and health care personnel as source control in accordance with CDC recommendations on Interim Infection Prevention and Control.

Disclaimer: Neither will it be intended for medical purposes, nor will it be considered as medical devices.

【Caution】 Do not use if the packaging is broken. The face mask is one time use only.

【Storage】 Store in a cool, dry place. Keep away from direct sunlight.

【Manufacturer Information】 EcoGuard, Inc. 700 S Battleground Ave#103, Grover, NC 28073

【Place of Origin】 United States



Produced in clean room



Hypoallergenic



Adjustable nose strip



Comfortable



Four layers protection



Soft and comfort

80mmx140mm